


Fumonisin

B1, B2, and B3

Use the Romer Labs® Fumonisin Clean-up columns for easy Fumonisin analyses.



3-step clean up

MultiSep®211 Fum




- SPE column
- Stability: 18 months
- Recoveries: >80%
- Repeatability: <5%
- Time for clean up: <20 minutes
- For corn, corn products and grains



Procedure for Fumonisin B1, B2, and B3 Clean up

3-step clean up

MultiSep® 211 Fum principle



Order Information:

MultiSep®211 Fum
Order No.:COCMU2211, Box of 25

Standard:
Fumonisin Mixture
(FB1 and FB2)
Order No.:002006